


Protein A resin selection guide



THIS IS AN INTERACTIVE DOCUMENT
Clicking on the product name will take you directly to the product specifications. If you require support or have any questions, please [talk to a specialist](#).



All Praesto Protein A resins are available in the following pre-packed/pre-qualified columns or bulk options.



OPUS® RoboColumn®

Volume:
0.5 ml - 0.6 ml

Use:
For use with robotic workstations for HTPD work.



OPUS® MiniChrom Column

Volume:
0.2 ml - 10 ml

Use:
For process development and parameter screening as well as small scale purification or sample preparation.



HT Column

Volume:
1 ml - 5 ml

Use:
For efficient resin screening for further optimization and verification.



OPUS® 5 - 80R Column

Volume:
0.5 L - 150 L

Use:
Designed to meet the chromatography requirements of larger 1000L and 2000L single-use bioreactors.



Bulk Resin

Volume:
10 ml - 10 L