

PRODUCT DECLARATION

Purolite® A430MR

This declaration was prepared by Purolite Corporation
in February, 2017.

Prepared by:
Doina Florea –
Quality & Regulatory Manager

Cosmin Catanas –
Quality Assurance Engineer



Purolite® Life Sciences

Victoria, Romania February 2017

PUROLITE S.R.L.

Str. Aleea Uzinei, Nr.11, 505700
Victoria, Jud. Braşov, România
T: +40 268 243 004 to 243 010
F: +40 268 243 002 and 243 003
office@purolite.ro
www.purolite.com
Nr. Reg. Com.: J08/446/1995
Cod Fiscal: RO 6039433

PRODUCT INFORMATION**General**

Product Name:	Purolite® A430MR
Manufacturer Code:	90612
Product Type:	Active Pharmaceutical Ingredient
Principal Applications:	Cholestyramine Resin
Manufacturer Site:	SC Purolite SRL
Chemical Name:	Poly(trimethylammoniomethyl styrene chloride- co-divinylbenzene) (x:y)
Ph. Eur. Name:	Cholestyramine
USP/NF Name:	Cholestyramine Resin
JP Name [JP 16]:	Cholestyramine
CAS Number:	11041-12-6

Regulatory status

- Drug Master File (DMF) registered on various markets: USA, Europe, Japan, Canada and Mexico.

Certifications / Approvals

- ISO 9001:2008 Certificate
- US-FDA Approval
- EU-GMP Certificate
- Kosher Certificate
- Halal Certificate

Material origin:

- Purolite® A430MR or any of its raw materials are of synthetic origin. It is manufactured using no raw materials or additives of human, animal and vegetal origin;
- Purolite® A430MR [Comprehensive Declaration](#) includes:
 - TSE/BSE free statement;
 - GMO free statement;
 - Allergens free statement;
 - Aflatoxin free statement;
 - Gluten free statement;
 - Residual solvents statement;
 - Residual metals statement;
 - Elemental impurities statement;
 - Genotoxic impurities statement;
 - Dioxins free statement;
 - Melamine and glycerol free statement.

Manufacturing

- Purolite® A430MR is manufactured through chemical synthesis;
- Purolite® A430MR is manufactured as per a systematic [Manufacturing Process](#);
- Purolite® A430MR batch size:

1500 kg packed in 50 kg drums;
1440 kg packed in 90 kg drums.

Packaging details


- Purolite® A430MR is packed in double polyethylene lined fibre drums; the internal liner is tie sealed and the outer one heat sealed. Fibre drum is sealed.

Pharmaceutical fibre drums		
DRUM TYPE	DIMENSIONS cm (in)	QUANTITY / NET WEIGHT
Medium	49 x 84.2 (19 x 33)	50 kg
Large	58.3 x 91.4 (23 x 36)	90 kg

Shelf-life and storage conditions

- When stored under proper storage conditions (in original packaging without extended exposure to temperature extremes or direct sunlight), Purolite® A430MR retain properties for a period of five (5) years from the date of manufacture;
- Purolite SRL has performed stability studies in order to establish the shelf life of the resin; stability data are available to support the assigned shelf life of five (5) years.

Manufacturing process flow chart
Cholestyramine – Purolite® A430MR



Date issued: **February 2017****DECLARATION**

Component Name: Purolite® A430MR, Cholestyramine
Manufacturer Catalogue No.: 90612
Component Grade: USP / Ph. Eur. / BP

Declaration statement for Purolite® A430MR

Purolite certifies that:

- Purolite® A430MR is a powdered, fully synthetic, polystyrenic, strong base anion exchange resin in chloride form, fully conforming to the USP monograph of Cholestyramine and to the Eu. Phar. / BP monograph of Cholestyramine;
- Purolite® A430MR is produced by our subsidiary PUROLITE S.R.L. in their GMP certified manufacturing plant located in Victoria, Romania, which is FDA inspected and compliant;
- Purolite Ltd and Purolite S.R.L. operate a Quality Management System according to BS EN ISO 9001:2008 that is certified by BSI.

We also certify that **Purolite® A430MR:**

- Is manufactured using no raw materials or additives of human, animal and vegetal origin, and it is consequently BSE and TSE safe and aflatoxin free;
- Is manufactured using no genetically modified organisms or products, therefore is GMO free;
- Is manufactured using NO raw materials or additives containing allergenic materials according to the ALBA list and the European Directive 2003/89/EC (gluten, crustaceans, eggs, fish, peanuts, soya beans, milk, nuts, celery, mustard and sesame seeds, and other products thereof);
- Is manufactured only with Methanol, a Class 2 Solvent; the residual methanol content of this product is less than 1,000 ppm, well below 3,000 ppm allowed;
- Is manufactured without the use of metal catalysts; as a result, material is free of residual metals deriving from catalysts;
- Is compliant to ICH Guidelines Q3D on elemental impurities, to USP/NF chapter <232>, current version;
- Is not manufactured with compounds with genotoxic properties, and genotoxic compounds are not created during the production process; as a result, the material is free of genotoxic impurities;
- Is not manufactured with dioxins, and dioxins are not created during the production process; as a result, the product is dioxin free;
- Is not manufactured with melamine and glycerol material and melamine and glycerol material is not created during the production process; as result, the product is melamine free and glycerol material free.

Doina FLOREA

Quality and Regulatory Manager

Americas

150 Monument Road
 Bala Cynwyd, PA
 19004
 T +1 800.343.1500
 T +1 610.668.9090
 F +1 484.384.2751
 Americas@purolite.com

Europe

Llantrisant Business Park
 Llantrisant
 Wales, UK
 CF72 8LF
 T +44 1443 229334
 F +44 1443 227073
 Europe@purolite.com

Asia Pacific

Room 707, C Section
 Huanglong Century Plaza
 No.3 Hangda Road
 Hangzhou, Zhejiang, China 310007
 T +86 571 876 31382
 F +86 571 876 31385
 AsiaPacific@purolite.com



Australia
 Bahrain
 Brazil
 Canada
 China
 Czech Republic
 France
 Germany
 India

Indonesia
 Israel
 Italy
 Japan
 Jordan
 Kazakhstan
 Korea
 Malaysia
 Mexico

Morocco
 Poland
 Romania
 Russia
 Singapore
 Slovak Republic
 South Africa
 Spain
 Taiwan

Tunisia
 Turkey
 UK
 Ukraine
 USA
 Uzbekistan



Purolite—the leading manufacturer of quality ion exchange, catalyst, adsorbent and specialty high-performance resins—is the only company that focuses 100% of its resources on the development and production of resin technology.

We're ready to solve your process challenges. For further information on Purolite® products and services, visit Purolite.com or contact your nearest Technical Sales Office.

©2017 Purolite
 All rights reserved
 PSI-01-A430MR-Rev.01/February 2017